

# Copper and Copper Alloys

**EXTEC®**      Δ      **SIMPLICITY**

## *Color Guide to Materials Preparation*

Copper and Copper alloys are difficult to prepare because of their softness and the secondary phase which typically is an iron oxide that can be dislodged (pulled out) and/or scratch the polished surface during preparation. The key to polishing is to reduce the amount of smeared deformation caused by grinding while retaining the inclusions during final polishing.

### **Hardware** I

1. Extec Labcut 250B Abrasive Cutting Machine ([www.extec.com/labcut250B](http://www.extec.com/labcut250B))
  2. Extec Labpress 40 Automatic Mounting Press ([www.extec.com/labpress40](http://www.extec.com/labpress40))
  3. Extec Labpol 12-3DI Auto Polisher/Grinder ([www.extec.com/labpol12-3DI](http://www.extec.com/labpol12-3DI))
- or
4. Extec Labpol 12 Auto Polisher/Grinder ([www.extec.com/labpol12](http://www.extec.com/labpol12))

### **Sectioning** II

A rubber – silicon carbide abrasive blade is preferred for medium hardness materials.

### **Mounting** III

Compression mounting with phenolics is typically used like our number 14505 black epoxy-mounting compounds.

## Grinding/Polishing Method -- Copper and Copper Alloys

Surface	Code	Abrasive/ Type Size	Lubricant	Code	Pressure (psi)	Time	Wheel Speed	Head Speed/ Direction
Coated Abrasive	VI	320 SiC	Water		5 psi	60 seconds	120 rpm	60rpm/Comp
Coated Abrasive	VI	400 SiC	Water		5 psi	60 seconds	120 rpm	60rpm/Contra
Coated Abrasive	VI	600 SiC	Water		5 psi	60 seconds	120 rpm	60rpm/Contra
Coated Abrasive	VII	P1200 SiC	Water		5 psi	60 seconds	120 rpm	60rpm/Comp
Coated Abrasive	VII	P2400 SiC	Water		5 psi	60 seconds	120 rpm	60rpm/Contra
S-Plan	VII	1um Diamond	Water Soluble Diamond Extender	I	5 psi	3 minutes	120 rpm	120rpm/Comp
S-Plan	VIII	0.25 um Diamond	Water Soluble Diamond Extender	I	5 psi	2 minutes	120 rpm	120rpm/Comp
Multi Cloth or Chemic Cloth	VIII	Colloidal Silica .06um	The last 10 seconds wash with Distilled Water		5 psi	90 seconds	120 rpm	120rpm/Contra

fm: 0420-20030325

